

Large-scale Electronic Balance



Introduction:

Large-scale electronic balance has fluent and novel appearance as well as new design, which bring faster weighing response. Unique metal base and mechanical structure guarantee higher stability and stronger anti-interference performance.

AC/DC exchangeable for portability and field operation.

Besides, we can also customize special specifications upon customer's request.

Feature:

- * The unique appearance design makes the product have high stability, and the backlit LCD screen makes the operator's vision more comfortable.
- * With AC and DC power supply, valve regulated lead-acid battery.
- * With a variety of measurement unit conversion, count, percentage, multi-point calibration.
- * The built-in RS232C output interface can be used with computers, printers and other equipment.

Technical Parameters:

Model	Capacity	Readability	Power Supply	Scale Size (mm)	External Size (W*D*H) mm	Net Weight	Package Size (W*D*H) mm	Gross Weight
BE40000L	0-40000g	1g	100~240V, 50/60Hz	500*400	530*300*680	11kg	870*465*225	14kg
BE50000L	0-50000g							
BE60000L	0-60000g							
BE100000L	0-100000g							
BE150000L	0-150000g							
BE200000L	0-200000g							
BE300000L	0-300000g							
BE50000L-10	0-50000g	10g						
BE60000L-10	0-60000g							
BE100000L-10	0-100000g							
BE150000L-10	0-150000g							
BE200000L-10	0-200000g							
BE300000L-10	0-300000g							
BE500000L-10	0-500000g	50g						
BE100000L-50	0-100000g							
BE150000L-50	0-150000g							
BE200000L-50	0-200000g							